HUMAN CARE
Your Life. Your Way.

FL 180/FL320

Maintenance Manual/
servicemanual

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Content:

Product Description	3
Maintenance, inspection and test	4
Maintenance checklist	5
Replacement of electrical parts	6

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Product Description

The intended use of the mobile lift is to transfer an individual from one resting surface to another (such as a bed to a wheelchair). Moving a person in a sling over ANY distance is not recommended.

The lift is intended for use in health care, intensive care and rehabilitation.

This maintenance manual is valid for FL180 and FL320.

Part numbers:

FL180:

70967-NA

70967-AU

70967-GB

70967-JP

70967-EU

FL320:

71967-NA

71967-AU

71967-GB

71967-EU

WARNING!

Service work contained in this manual should only be undertaken by Human Care trained and authorised service agents. Only Human Care supplied parts should be used. Failure to comply with the above two items will void warranty and could result in further unnecessary damage to the product and possible extreme safety hazard to the patient or caregiver.

PURPOSE

This document will be made available to Human Care staff, distributors, agents and authorised Human Care service providers. It is not for publication or distribution to the public or clients who may use the information to allow other companies or individuals to repair or modify Human Care products.

This document and its contents are intended to be used in the performance of service on Human Care produced products, i.e. FloorLine LTC range series products. This document also serves as a reference for office and sales staff in the performance of their duties, i.e. collecting service call information to facilitate the expedient diagnosis of any problems and problem solutions.

This document should be read in conjunction with Human Care's Operating and Assembly Guides, Instruction Manuals and Spare Parts Catalogues. This document is copyright to Human Care and must not be reproduced in any form without written permission of Human Care.

MAINTENANCE, INSPECTION AND TEST

AS a guideline it is recommended to perform an annual inspection of the lift and its accessories. The examination and test should be conducted according to the recommendations and procedures below. If any damage or malfunction is suspected the lift has to be redrawn from operation and Human Care should be contacted.

Inspection list

1. **SPREADER BAR:** Check the spreader bar for freedom of rotation and swing. Check for wear on the central pivot. Check for firm attachment to the boom.
2. **BOOM:** Check the attachment of the boom to the mast. Make sure there is only minimal side movement of the boom and the boom is free to rotate on the boom mounting.
3. **MAST:** Check the operation of the mast-locking device. Make sure the mast fully engages into the socket. Check the bottom actuator unit mounting.
4. **ACTUATOR:** The actuator should require no maintenance other than checking for correct operation and listening for unusual noise.
5. **POWER PACK:** Check the function of emergency stop button and emergency down/up.
6. **BATTERIES:** The batteries are housed in the power pack and should not require maintenance other than regular charging as detailed in the charging instructions. Check that the connections remain clean.
7. **LEG ADJUSTMENT:** Check the legs operate in both full extensions (inward/outward).
8. **CASTORS:** Check all castors for firm attachment to the legs. Check for free rotation of the castor and the wheels.
9. **CLEANING:** Clean with ordinary soap and water and/or any hard surface disinfectant. Harsh chemical cleaners or abrasives should be avoided as these may damage the surface finish of the lift. Avoid wetting any of the electric parts.
10. **BASE AND WHEELS:** Ensure base is even and level (all four wheels are on the floor).
11. **SLINGS:** Check for wear and fraying.
12. **LUBRICATION:** Oil pivot joints, including mast and boom connections, pedal assembly, and spreader bar joint (only if required).

Maintenance checklist

ENG

Check to be performed		OK	Not OK	Description
Hanger bar	Damage, scratches, function			
Boom	Damage, scratches, function			
Mast	Damage, function			
Spreader Actuator	Damage, noise, function			
Lifting actuator	Damage, noise, function			
Castors	Damage, locking function, attachment			
Slings	Wear, damage			
Emergency stop	Function			
Emergency lowering mechanical	Function			
Emergency lowering control box	Function			
Battery pack	Damage			
Base and castors	Even level			
Hand control	Damage			
Hand control cable	Damage, routing			
Cabling actuators-control box	Damage, routing			

REPLACEMENT PARTS

In order to maintain operational safety and the right to claim under warranty, only original Human Care replacement parts may be used! Please quote part numbers for each item as per the Human Care Spare Parts Catalogue.

ELECTRICAL COMPONENTS

The corresponding replacement parts can be obtained from Human Care. When ordering please specify – TYPE & ITEM. The necessary details are found on the type plate (as shown below) on each electrical component.

Safety information on replacing electrical components

This is a Class 2 electrical product.

CAUTION! Before commencing any work on electrical equipment, always unplug the mains cable from the electrical socket.

Any work and/or repairs to the electrical equipment may only be carried out by service engineers, qualified technicians or authorised electricians in compliance with all the relevant country of use safety standards.

All electrical drive system components (actuators, handset, and control box) are maintenance-free and must not be opened.

In the event of a malfunction, the corresponding component must always be replaced in full.

Change Spreader Actuator FL180

- Disconnect cables between control box and actuators
- Remove the mast from the base

- Place the base piece of the lift upside down.
- Remove the plastic caps from the nuts and loose the screws holding the spreader actuator
- Connect the actuator power cable (p.n 89610) to the new LA23 actuator (p.n. 70973). Lock the cable with the cable lock (p.n 89611).
- Assemble the LA23 to the metal tubes. Cover the nuts and bolt with plastic caps. Place the washers in the order shown in the picture below.

- Put the mast back into the base and connect all cables
- Check the function of the new actuator.

Change Spreader Actuator FL320

- Disconnect cable between control box and actuator
- Remove the mast from the base

- Place the base piece of the lift upside down
- Open the plastic cover and loose the screws holding the spreader actuator.
- Assemble the LA28 p.n 71230 to the metal tubes.

- Assemble the plastic cover
- Put the mast back into the base and connect all cables
- Check the function of the new actuator

Change Lifting Actuator/Control box

- The lifting actuator and the control box have to be exchanged at the same time if calibration cannot be performed onsite.
- Make sure that the lifting actuator and the control box have been paired at Human Care before exchanging
- Remove cabling between control box and actuators
- Remove hand control
- Remove battery pack
- Remove the control box (loosen screw)
- Attach the new control box the mounting plate
- Use the screw that comes with the mounting plate to attach the CBJ-C control unit (part number FL180: 70971, FL320: 71232) to the mounting plate.

- Connect the BAJ1 battery pack (part no. FL180: 70975, Part no. FL320: 71234).

- Connect cable to the actuator
- Attach the hanger bar
- Verify the lifting function of the actuator.

- Remove hanger bar.
- Remove the plastic caps from the nuts holding the lifting actuator and loose the screws.
- Remove the LA34 (FL180)/LA44 (FL320) actuator.
- Attach the LA34 (FL180)/LA44 (FL320) actuator (part no. FL180: 70972, Part no. FL320: 71228) to the tower and the lift arm. Cover the nuts and bolt with the plastic caps. Place the washers in the order shown in the picture below.

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