
## **Activity Sling Fitting Protocol**



Client Informat	Name:		Room:	
Facility:		Lift Used:	Hanger Bar Use	ed:
Sling Model:	Name/Article Number:	Fabric	Size	
	Loops for back:	Green Strap Yellow Strap	Red Straps	Extension Loops Back: 21180-10 (10cm/3.9") 21180-15 (15cm/5.9") 21180-20 (20cm/7.8") 21180-25 (25cm/9.8") 21180-35 (35cm/13.7")

\*Please note that the leg placement information box below is not relevant for the Activity Sling. \*Only 2-point suspension and only arm placement on the outside of the sling is possible for the Activity Sling.



Notes:

## Size guidelines:

KG		↓ <u>C</u>			JS JM JL S M L XL XL XXX
User weight kg	User length cm	User back height* cm	User hip width* cm	Recommended size*	Färg
15-20	110-130	69-79	25-30	Junior Small (JS)	JS
20-35	120-150	75-80	29-32	Junior Medium (JM)	M
35-50	130-160	81-91	31-36	Junior Large (JL)	L
45-60	140-170	87-97	35-38	Small (S)	S

 $^{*}$  To measure the client back height, measure from the clients coccyx to the top of the head.

\* To measure the client hip width, have the client sitting down and measure in a straight line (not the circumference).

\* These sizes are recommendations only. A correct fitting must always be carried out by clinically knowledgeable personnel using the sling's user manual. Please be aware that some slings are not available in all sizes. You can find drawings on each sling in the user manual, delivered with the product and available on our websites.

## Fitting Protocol Completed by:

i

Name:	
Date:	Position/Title:

Important! This fitting protocol does not replace the user manual. Read all instructions before use!

www.humancaregroup.com