



Cable-free

Wide backstrap

Three colours



Our bestselling classic

neXus 3, is our bestselling rollator with a wide backstrap and robust wheels, making it a perfect everyday-life-companion.



Easy to fold

The best selling neXus series rollators have innovated the way we think about walking aids. With the first-ever patented cross-folding frame, it set the standard for future walking aids design. The neXus 3 rollator features a cross-folding frame, making it exceptionally easy to fold. The soft zippered bag folds together while remaining on the rollator. It can stand independently in a folded position with the bag on the rollator. The cross-folding design also enables the rollator to be partially folded to pass through narrow spaces.

braking system also ensures years of trouble-free use. We have had anatomy and strong visual communication in mind when designing the brake handles. They are comfortable, and the red ball makes the brake handles visible. The cross-folding structure enables superior walking ergonomics with ample walking space. The padded, contoured, and cushioned seat offers a comfortable resting spot. The backstrap can be flipped to allow for two-position seating.

Available in black, blue and red

The neXus 3 is available in three different sizes; neXus 3 Super Low, neXus 3 Low, and neXus 3 Standard. All sizes are available in Black, Blue and Red (Burgundy) colour.




Cable-free and wide backstrap

The cable-free and anatomic brake handles increases safety and comfortability. No cables catch on objects around the home. The cable-free



Features




Padded Backrest

Back comfort while seated. Can be flipped both ways and is included with the rollator.



Foam Cushioned Seat

Offers a comfortable resting spot while walking.



Anatomic Brake Handles

The anatomic brake handles are easy and comfortable to use.



Tray

Practical for resting lightweight items. Maximum load 2 kg/4.4lb. Part no. 93316.



Hard Basket

Basket ideal for shopping. Maximum load 4.5kg/10 lb. Part no. 4510.



Cane or Umbrella Holder

Simplifies transport of a cane or umbrella. Part no. 93317.



Cup Holder

Attaches in seconds without tools. Fits standard size bottle and cups. Part no. 90394.



Curb Climper

Allows the rollator to be slightly lifted onto a curb. Part no. 93315.



Identification Tag Green

Attaches to the rollator in seconds. Part no. 90393.



Identification Tag Canada

Attaches to the rollator in seconds. Part no. 90392.



Soft Zippered Bag

Soft bag that folds with the mobility aid. Maximum capacity 5kg/11lb. Part no. 93280.

Technical specifications

Part Name	Part no.	Handle Height (Min - Max)	Seat Height	Overall Width	Width between Handles	Distance between Rear Wheels	Wheel Size	Product Weight
neXus 3 Super Low		78 cm - 86 cm	48 cm	58 cm	45 cm	47 cm	15 cm x 3 cm	7.6 kg
Black	470018	30.5" - 34.1"	19"	23"	17.5"	18.5"	6" x 1.2"	16.9lb
Blue	470018B							
Red (Burgundy)	470018R							
neXus 3 Low		78 cm - 93 cm	56 cm	58 cm	45 cm	47 cm	15 cm x 3 cm	7.7 kg
Black	470022	30.5" - 36.5"	22"	23"	17.5"	18.5"	6" x 1.2"	17.1lb
Blue	470022B							
Red (Burgundy)	470022R							
neXus 3 Standard		83 cm - 101 cm	61 cm	58 cm	45 cm	47 cm	15 cm x 3 cm	7.8 kg
Black	470024	32.5" - 39.8"	24"	23"	17.5"	18.5"	6" x 1.2"	17.3lb
Blue	470024B							
Red (Burgundy)	470024R							

Maximum User Weight: 200 kg / 441 lb



The neXus 3 is available in three different colours; Black, Blue, and Red [Burgundy].



*We improve freedom of movement
and quality of care for people.*



Human Care Canada
10-155 Colonnade Road
Ottawa, ON K2E 7K1, CANADA
+1 613 723 6734
info.ca@humancaregroup.com

www.humancaregroup.com