



Alpha Stand-up Support Magic & Standard

Always
ready
to support




With the Alpha Stand-up Supports you can securely hold on to, and guide, the client when standing up, sitting down and particularly during transfers. When being used, the simple handling can easily be combined with a back-friendly work method and elements from nursing concepts. The Alpha Stand-up Supports can be combined with other aids in various ways, depending on the clients needs and functional requirements. The extremely stable carrying handles provide additional safety.



The Alpha Looped Belt, PM-3010, can be used optionally for clients with limited body tension. It enables also especially the safe transfer for clients that sit in very narrow [wheel]chairs where the flexible support blet of the Stand-up Support Magic cannot be applied to the client from behind.

Technical specification

| Product Name | Art. nr. | Dimensions | Material | Max lifting weight |
|-------------------------|----------|----------------------------|-----------------------|--------------------|
| Alpha Stand-up Standard | PM-3012 | 60 x 20 cm 23,6" x 7,8" | Polyamide, latex free | 200 kg / 440 lb |
| Alpha Stand-up Magic | PM-3015 | 60 x 20 cm 23,6" x 7,8" | Polyamide, latex free | 200 kg / 440 lb |



These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.