


# Alpha Stand-up Aid Pegasus

Compact  
Manual  
Light



The Alpha Stand-up Aid Pegasus is used to support a person after a fall to assist them from a supine position from the floor to a comfortable seated position. Depending on the clinical condition, this can be performed by themselves or with help of a care giver. It is very compact, safe and lightweight to handle.

This raising aid is ready for immediate use at home or as an ideal aid for utilization to provide support and help. The seat of the Stand-up Aid Pegasus has a minimum receiving height of 8 cm. The raising aid reaches a seat height of 42 cm.

## Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Stand-up Aid Pegasus	PM-4222	Standing 36 x 64 cm included back 50 cm  /	Drive mechanism: - high grade stainless steel and spring steel Seat: Polypropylene Back: Polypropylene Handle: High grade stainless steel 1.4301 Handle box: Plastic white Frame: High grade stainless steel 1.4301	135 kg / 298 lb

   = 135 kg / 298 lb

**PETERMANN**  
part of the Human Care Group

These products are distributed by Human Care AB for Petermann GmbH.



Visit [www.humancaregroup.com](http://www.humancaregroup.com) for the latest product information and contact details to national offices and partners.