

fellow classic



The Fellow Classic is a high quality rollator ideal for both indoor and outdoor use. Safe and reliable, it is available in five different models - standard, low, high, wide and narrow. The soft braking system is ideal for those users with limited hand mobility.

Unique Design

The Fellow Classic features a unique fork construction that makes it extremely easy to maneuver, especially in smaller spaces. The easy-to-use folding system makes storage and transport easy and the solid wheels prevent build up of dirt and rust allowing them to operate and move more smoothly.

Ergonomic Grips

The Fellow Classic's features handles that are ergonomically shaped to distribute the weight put on them evenly. Combined with the soft braking system, this ensures that minimal hand pressure is required to lock the wheels in place. The forward facing handles allow the user to walk in a more upright

position which is easy on both the back and shoulders.

Practical Accessories

The Fellow Classic comes with an assortment of accessories to make day-to-day life easier. From a tray to a cane holder, you can easily find an accessory to suit your needs. The grid easily removes for cleaning and maintenance.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone: +46 8 665 35 00
Fax: +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone: +61 1300 132 736
Fax: +61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features




Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 08806



Soft Bag

Soft carry case ideal for personal items. Part no. 08808



Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



Mesh basket

Ideal for those great shopping adventures. Part no.: 88051 (low), 88052 (narrow, wide, high), 08805 (standard)



Light kit

Illuminate your way while walking. Part no. 91091



Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603



IV Stand

Part no. 08814



One-Hand Brake

Allows control of brakes with one hand. Must be ordered with rollator. Part no. 10677 (right), 10678 (left)



Slow-Down Brake

Prevents the rollator from moving too fast. Part no. 10625 (set), 88121 (right), 88122 (left)



Backrest

Backrest for extra comfort while seated. Part no. 10668



Curb Climber

Helps raise the rollator onto a curb. Part no. 10661

Technical Specifications

	Standard	Low	Wide	Narrow	High
Part no.	88009	88109	88309	88509	88609
Handle Height (Min - Max)	77cm - 95cm / 30.3" - 37.4"	68cm - 81cm / 26.8" - 31.7"	77cm - 95cm / 30.3" - 37.4"	77cm - 95cm / 30.3" - 37.4"	82cm - 100cm / 32.3" - 39.4"
Seat Height	59cm / 23.2"	51cm / 19.9"	60cm / 23.6"	60cm / 23.6"	64cm / 25.2"
Overall Width	63cm / 24.8"	62cm / 24.4"	69cm / 27.2"	60cm / 23.4"	63cm / 24.8"
Overall Length	69cm / 27.2"	60cm / 23.6"	69cm / 27.2"	69cm / 27.2"	69cm / 27.2"
Width between Handles	44cm / 17.3"	44cm / 17.3"	50cm / 19.7"	44cm / 17.3"	44cm / 17.3"
Wheel Size	20cm x 4cm / 7.7" x 1.7"	20cm x 4cm / 7.7" x 1.7"	20cm x 4cm / 7.7" x 1.7"	20cm x 4cm / 7.7" x 1.7"	20cm x 4cm / 7.7" x 1.7"
Product Weight	8kg / 17.6lb	7.7kg / 16.9lb	8.7kg / 19.2lb	8.1kg / 17.8lb	8.1kg / 17.8lb
Maximum User Weight	150kg / 330.7lb	130kg / 286.6lb	150kg / 330.7lb	130kg / 286.6lb	130kg / 286.6lb

This product is manufactured in accordance with all relevant ISO standards.



Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup



youtube.com/channel/UC7fjN3zTjRxOP1IDf2Da0

humancaregroup.com