



Lightweight
Easy to Fold
Compact



The Parade is the perfect choice when looking for a lightweight, versatile rollator. Its design ensures easy manoeuvrability even in small spaces. Despite its lightweight frame, it is still very stable and offers outstanding value for money.



Easy-to-Fold Frame

The Parade features an integrated handle in the seat that makes folding easy. Simply slide the seat over and pull on the handle. It's that easy! The Parade stands independently, even in the folded position. It is a lightweight rollator that is easy to manoeuvre, even in smaller spaces. The Parade has a maximum user weight of 140kg/309lb.

Practical Accessories

The Parade comes with an assortment of accessories to make day-to-day life easier including baskets and backstraps. It is the perfect option for those outdoor adventures or for just moving in and around the home.



Ergonomic Grips

The Parade features ergonomic handles designed to provide users with reduced hand functions the ability to easily lock the wheels in place and comfortably maneuver the rollator while walking.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features



Tray

Practical for resting lightweight items.
Maximum load 2kg/4.4lb.
Part no. 10664




Curb Climber

Helps raise the rollator onto a curb.
Part no. 10661



Mesh basket

Ideal for those great shopping adventures.
Maximum capacity 4.5kg/10lb.
Part no.: 10662



Light kit

Illuminate your way while walking. Part no. 91091



Cane holder

Simplifies transport of a cane or umbrella.
Part no. 10675



Reflector Kit

Includes 30 self-adhesive reflectors.
Part no. 91090



Backrest

Supports the back while resting. Part no. 10621



Security Lock

Flex cord that protects rollator from theft while unattended.
Part no. 91092

Technical Specifications

Low	
Part no.	85500
Handle Height (Min - Max)	74cm - 95cm / 29.1" - 37.4"
Seat Height	60cm / 23.4"
Overall Width	59cm / 23.2"
Overall Length	69.5cm / 27.4"
Width between Handles	43cm / 16.9"
Distance between Rear Wheels	50.8cm / 20"
Wheel Size	20cm x 4cm / 7.7" x 1.7"
Product Weight	8.8kg / 19.4lb
Maximum User Weight	140kg / 309lb

This product is manufactured in accordance with all relevant ISO standards.



Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup



youtube.com/channel/UC7fjiNh3zTjRxOP1IDf2Da0

humancaregroup.com