



The Staffan RA is a highly maneuverable rollator primarily for indoor use. It is only 59cm (23.2") wide, making it an ideal option for moving around tighter spaces in and around the home.

Compact Design

The Staffan RA features fixed rear wheels and brakes for better stability and security as well as a narrow overall width that makes it ideal for use in the home. Its easy-to-fold frame makes transport and storage easy. Designed for users that have to deal with walking while suffering from arthritic or neurological conditions, the Staffan RA will allow for maximum stability while performing everyday functions.

Unique Handles

The Staffan RA features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with

integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight. The Staffan RA can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Adjustable Forearm Supports

The Rebel RA features very unique forearm supports that are adjustable in both height and direction. This ensures that the forearm supports can be positioned in the best possible way for user comfort and stability.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features




Reflector Kit

Includes 30 self-adhesive reflectors.
Part no. 91090



Security Lock

Flex cord that protects rollator from theft while unattended.
Part no. 91092




Basket

Included with rollator in Canada. Ideal for shopping Part no. 40601




Light kit

Illuminate your way while walking. Part no. 91091




Cane holder

Simplifies transport of a cane or umbrella.
Part no. 40603



One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Part no. 10653 (low), 10654 (high)



Oxygen Holder

Holds oxygen tube firmly in place. Part no. 40615

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.



Technical Specifications

	Staffan RA Low 59	Staffan RA High 59
Part no.	42221	42231
Handle Height (Min - Max)	86cm - 108cm / 33.8" - 42.5"	107cm - 128cm / 42.1" - 50.4"
Overall Width	59cm / 23.2"	59cm / 23.2"
Wheel Size	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"
Product Weight	10.3kg / 22.7lb	10.6kg / 23.4lb
Maximum User Weight	90kg / 198.4lb	125kg / 275.6lb