Carl-Oskar TW












The Carl-Oskar Trolley features a wide, comfortable forearm support platform table that allows pain sensitive users to comfortably rest their forearms while the front-facing brake handles assist to maintain better posture.



Outstanding Stability

The Carl-Oskar TW features a reactive frame that remains stable over virtually any surface. With extra large ball-bearing mounted wheels and a comfortable and firm table, the Carl-Oskar TW is the ideal option for pain sensitive users that find it difficult to use regular trolley walkers.






The Carl-Oskar TW features unique frontfacing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight.

Optional Accessories

The Carl-Oskar TW offers a vast assortment of accessories including oxygen tank holders that allows the walker to function for a number of users. The Carl-Oskar TW easily folds once the table is removed.



HUMAN CARE SWEDEN (HQ) HUMAN CARE CANADA

Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613 723 6734 Fax: +613.723.1058 info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

Accessories/Features



Tray
Practical for resting
lightweight items.
Maximum load 2kg/4.4lb.
Part no. 10601



Splash GuardProtects legs and clothes.
Part no. 10608



Reflector Kit Includes 30 selfadhesive reflectors. Part no. 91090



Slow-Down Brake Prevents the rollator fro moving too fast and enables a controlled walking speed for the user. Part no.: 10624



Light kitIlluminate your way while walking. Part no. 91091



Security Lock
Flex cord that protects
rollator from theft while
unattended.
Part no. 91092



Cane holder
Simplifies transport of a cane or umbrella.
Part no. 10603



One-Hand Brake
Allows user to control
hand functions with one
hand (installed left or
right). Special order item.
Part no. 10653



Curb Climber Helps raise rollator onto a curb. Part no. 10612



Oxygen Holder Holds oxygen tube firmly in place. Part no. 10615



Basket Ideal for shopping. Part no. 10621

Technical Specifications

•				
	Carl-Oskar TW Low 67	Carl-Oskar TW Low 72	Carl-Oskar TW Nor 67	Carl-Oskar TW Nor 72
Part no.	13315	13415	13325	13425
Handle Height (Min - Max)	95cm - 113cm/ 37.4" - 44.5"	95cm - 113cm/ 37.4" - 44.5"	107cm - 126cm/ 42.1" - 49.6"	107cm - 126cm/ 42.1" - 49.6"
Overall Width	67cm/26.4"	72cm/28.3"	67cm/26.4"	72cm/28.3"
Width between Handles	47cm/18.5"	52cm/20.5"	47cm/18.5"	52cm/20.5"
Distance between Rear Wheels	50.8cm/20"	52.1cm/20.5"	50.8cm/20"	52.1cm/20.5"
Wheel Size	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"	30cm x 5cm/ 11.8" x 1.8"
Product Weight	14.5kg/31.9lb	14.6kg/32.2lb	14.7kg/32.4lb	14.8kg/32.6lb
Maximum User Weight	125kg/275.6lb	125kg/275.6lb	125kg/275.6lb	125kg/275.6lb

Note: Not all accessories are available for sale in all countries.