



The Rebel Trolley is a safe and reliable rollator with a firm table pad for extra support. Designed for users requiring that extra support, or those who find it difficult to maneuver an ordinary walker, the Rebel TW is the ideal rollator for both indoor or outdoor use.

Outstanding Stability

The Rebel TW is the ideal rollator for pain sensitive users suffering from arthritis, Parkinson's or many other neurological disorders. It features a wide, firm table support and large ball-bearing mounted wheels that ensure easy maneuverability and comfort. The Rebel TW can be used both indoors and outdoors and features a wide assortment of accessories for everyday use.

Unique Handles

The Rebel TW features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle

tube. The ergonomic design of the grips ensures even distribution of weight. The Rebel TW can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Reactive Frame

The Rebel TW has the same great frame as the regular Rebel which means that it remains stable even on uneven surfaces. This is especially helpful in rural areas where the ground is not even. The Rebel TW comes in two seat heights to accommodate most users.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street
Port Melbourne, VIC 3207
Phone : +61 1300 132 736
Fax : + 61 (03) 9645 4954
info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features



Tray
Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10601



Basket
Maximum load 4.5kg/10lb. Ideal for shopping. Part no. 10621



Reflector Kit
Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake
Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



Light kit
Illuminate your way while walking. Part no. 91091



Security Lock
Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder
Simplifies transport of a cane or umbrella. Part no. 10603



One-Hand Brake
Allows user to control hand functions with one hand (install left/right). Special order item. Part no. 10653



Oxygen Holder
Holds oxygen tube firmly in place. Part no. 10615



Curb Climber
Helps raise rollator onto a curb. Part no. 10612

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.



Technical Specifications

Description	Part No	Width between Handles	Overall Width	Handle Height (Min-Max)	Wheel Size	Product Weight	Maximum User Weight
Rebel TW 67 Low	23311	40.6cm/16"	67cm/26.4"	86cm-104cm / 33.8" - 40.9"	20cm x 4cm / 7.9" x 1.6"	12.8kg/28.2lb	125kg / 275.6lb
Rebel RA 67 Nor	23321	40.6cm/16"	67cm/26.4"	104cm - 122cm / 40.9" - 48"	20cm x 4cm / 7.9" x 1.6"	13.0kg/25.3lb	125kg / 275.6lb
Rebel RA 72 Low	23411	50.8cm/20"	72cm/28.3"	86cm-104cm / 33.8" - 40.9"	20cm x 4cm / 7.9" x 1.6"	12.9kg/28.4lb	125kg / 275.6lb
Rebel RA 72 Nor	23421	50.8cm/20"	72cm/28.3"	104cm - 122cm / 40.9" - 48"	20cm x 4cm / 7.9" x 1.6"	13.1kg/28.9lb	125kg / 275.6lb

Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup



youtube.com/channel/UC7fjInh3zTjRx0P1IDf2Da0