



The Carl-Oskar is the perfect rollator for the active user. Its durable frame is easy to maneuver and features ergonomic handles for a comfortable grip. With extra large wheels, it is the perfect rollator for those outdoor activities, even during snowy winters.

Outstanding Stability

The Carl-Oskar is one of Human Care's best known and best-selling rollators worldwide. It is safe and reliable and features extra-large 30cm/11" wheels that allow the rollator to roll very easily yet remain stable on virtually any surface. The Carl-Oskar is the ideal option for rural areas, even on snowy surfaces.

Unique Handles

The Carl-Oskar features unique angle adjustable handles that allow users to find the optimal handle position for comfort and stability. The ergonomic design ensures even distribution of weight for a stable grip. The Carl-Oskar can also be equipped with a one-hand brake to allow a user to control hand

functions on either the left or right hand side.

Reinforced Option

The Carl-Oskar is available in different overall widths to accommodate a wide variety of users. The reinforced models for bariatric users feature a weight capacity of 275kg/606lb. All models come with a full assortment of accessories to meet most needs and can easily be folded for storage or transport.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com


Accessories/Features



Tray
Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10651



Splash Guard
Protects legs and clothes. Part no. 10608



Reflector Kit
Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake
Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



Light kit
Illuminate your way while walking. Part no. 91091



Security Lock
Flex cord that protects rollator from theft while unattended. Part no. 91092




Cane holder
Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake
Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10652




Bag Kit
Includes one small (5kg capacity) and one large bag (10kg capacity). Part no. 91095 (for 62cm), 91096



Oxygen Holder
Holds oxygen tube firmly in place. Part no. 10615



Seat Pad
State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093



Backstrap and Basket
Included with rollator in Canada. Backrest Part no. 10621, 10655 (for 72cm) Basket Part no. 10605

Technical Specifications

Description	Part No	Width between Handles	Overall Width	Distance btwn Rear Wheels	Product Weight
Carl-Oskar 57	11150 11151 (Air Tires)	37cm / 14.6"	57cm / 22.4"	50.8cm / 20"	11.1kg / 24.5lb
Carl-Oskar 62	11250 11251 (Air Tires) 61225 (Reinforced)	42cm / 16.5"	62cm / 24.4"	50.8cm / 20"	11.3kg / 24.9lb 11.3kg / 24.9lb 13.2kg / 29.1lb
Carl-Oskar 67	11350 11351 (Air Tires) 61325 (Reinforced)	47cm / 18.5"	67cm / 26.4"	52.1cm / 20.5"	11.4kg / 25.1lb 11.4kg / 25.1lb 13.3kg / 29.3lb
Carl-Oskar 72	11450 11451 (Air Tires) 61425 (Reinforced)	52cm / 20.5"	72cm / 28.3"	52.1cm / 20.5"	11.5kg / 25.3lb 11.5kg / 25.3lb 13.4kg / 29.5lb

Description (All Models)	
Handle Height (Min - Max)	77cm - 95cm / 30.3" - 37.4"
Seat Height	64cm / 25.2"
Wheel Size	30cm x 5cm / 11.8" x 1.8"
Seat Size	31.8cm x 17.8cm / 12.5" x 7"
Depth folded	26cm / 10.2"
Maximum User Weight	Regular 150kg / 330.7lb Reinforced 275kg / 606.3lb

Note: Not all accessories are available for sale in all countries.

This product is manufactured in accordance with all relevant ISO standards.



Follow Us on Social Media



facebook.com/humancareworld



instagram.com/humancaregrip



linkedin.com/company/human-care-ab



pinterest.com/humancaregrip



youtube.com/channel/UC7fjInh3zTjRx0P1IDf2Da0