## staffan TW







Highly maneuvrable, this trolley walker is designed for indoor use only. Its narrow overall width makes moving around the home easy and its firm table provides comfort and support for those pain sensitive users.



### **Outstanding Stability**

The Staffan TW is the ideal rollator for pain sensitive users suffering from arthritis, Parkinson's or many other neurological disorders. It features a wide, firm table support and a narrow overall width that makes it easy to move in and around the home.



The Staffan TW features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight. The

Staffan TW can also be equipped with a onehand brake to allow a user to control hand functions on either the left or right hand side.

### Fixed Rear Wheels

The Staffan TW features fixed rear wheels and front swivel wheels for enhanced maneuvrability as well as a hand operated braking system for enhanced safety. Available in two heights, the Staffan folds easily for storage and travel.



#### **HUMAN CARE SWEDEN (HQ)**

Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

#### **HUMAN CARE AUSTRALIA**

Unit 4, 49 Bertie Street Port Melbourne, VIC 3207 Phone: +61 1300 132 736 Fax: +61 [03] 9645 4954 info.au@humancaregroup.com

## HUMAN CARE CANADA

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com

## HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

#### **HUMAN CARE UNITED STATES**

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

## **Accessories/Features**



Ideal for shopping. Part no. 40601



Reflector Kit Includes 30 selfadhesive reflectors. Part no. 91090



Cane holder Simplifies transport of a cane or umbrella. Part no. 40603



Light kit Illuminate your way while walking. Part no. 91091



Oxygen Holder Holds oxygen tube firmly in place. Part no. 40605




Security Lock Flex cord that protects rollator from theft while unattended. Part no. 91092



One-Hand Brake Allows user to control hand functions with one hand (installed left or right). Part no. 10653 (low), 10654 (high)

Note: Not all accessories are available for sale in all countries.





# **Technical Specifications**

	Staffan TW Low	Staffan TW High
Part no.	43211	43221
Handle Height (Min - Max)	86cm - 108cm / 33.8" - 42.5"	107cm - 128cm / 42.1" - 50.4"
Overall Width	62cm/24.4"	62cm / 24.4"
Wheel Size	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"
Product Weight	11.7kg/25.8lb	11.9kg/26.2lb
Maximum User Weight	90kg/198.4lb	125kg/275.6lb

Follow Us on Social Media









