



Compact
Maximum user weight
160 kg
Foldable



Ideal lift for many situations

The Diana II is an extremely stable, electrically operated lift with a compact, foldable frame so that it can be used in spaces with limited room to maneuver. Our extensive range of slings provides the required support.



The Diana II provides clients with the support they need when they are lifted from their bed, in/out of a chair or wheelchair, on/off the toilet and even from the floor. Thanks to its compact size the Diana II is ideal for use in small spaces, for instance, at home where there may be narrow doorways and in toilets/bathrooms. It is also very useful in healthcare institutions where other medical devices that are necessary can often clutter an area. The large, high-quality and durable twin casters ensure that the lift is easy to maneuver.

Various lifting options

The Diana II comes with a 2-point hanger bar or a manual 4-point spreader bar. The hanger bar is fitted with a loop attachment to which the loops of our slings can be securely attached. The manual 4-point spreader bar is fitted with mushroom shaped hooks to which the E-clips of our slings can be securely attached. The 4-point hanger can be used for positioning the client in three different client positions: recumbent, semi-recumbent and sitting. Our slings are attached using the unique E-clips.

The universal E-clip is a proven system for fastening a sling to a lift, thereby significantly reducing the risk of errors. The clip is placed around the mushroom shaped hooks of the spreader bar and closes with an audible 'click'. Our slings provide optimal support for extra secure, comfortable and safe transfers.



The leg spreading allows the lift to be placed around a chair/toilet to get as close as possible to a client for a more comfortable transfer. You can choose between operating the legs manually using foot pedals or by means of a fully integrated electrical control.



Mobility Aids



Healthcare Beds



Lifting Solutions



Hygiene



Convertible Chairs

Accessories



Full Body Sling E-clip Soft
Provides full support, ensuring the comfort and safety for every user and every situation. Attachment via secure E-clips or loops.
Article number 25010C-size



Hanger bar 4MSB
Manually operated 4-point hanger bar, with attachments for E-clip slings.
Article number 90558



Hanger Bar 2P
2-point hanger bar for use with a strap loop sling.
90567 (35 cm)
90557 (45 cm)
90568 (60 cm)



Soft padding for the hanger bar
Option to decide where to place the padding on the hanger bar.
Article number 50866

Use



Lifting in/out of bed
Provides support to clients during lifting from the bed and can help with positioning on the bed.



From and to a chair
Suitable for transferring to and from a chair or wheelchair and positioning the client in a chair.



Lifting from the floor
Provides support to clients while being lifted from the floor.



Toileting
Suitable for transfer to and from the toilet.

Features



Casters
2 twin castors, with brakes (Ø 100 mm, rear) and 2 twin castors, without brakes (Ø 75 mm, front).



Hand control
User-friendly hand control, easy-to-use, ergonomic design (1 x Manual Leg Adjustment MLA/1 x Electric Leg Adjustment ELA).



Foot pedal
Foot pedal for manual leg adjustment (MLA)



Push handle
Large push handle, which allows for grip at an ergonomic height, loop shaped so that clothing or slings do not get caught.

Technical specifications

Specifications	Diana II ELA, 90556 with Hanger Bar 2P	Diana II ELA, 90556 with Hanger Bar 4MSB
Length of the chassis	1280 mm	1280 mm
Width of the chassis	620 mm	620 mm
Height of the mast	1360 mm	1360 mm
Minimum width of the chassis (internal - external)	490 – 610 mm	490 – 610 mm
Maximum width of the chassis (internal - external)	1000 – 1120 mm	1000 – 1120 mm
Lifting range	720 – 1940 mm	720 – 1940 mm
Turning radius	1400 mm	1400 mm
Height of the chassis	95 mm	95 mm
Diameter of the casters	100 mm (rear) / 75 mm (front)	100 mm (rear) / 75 mm (front)
Total weight	45 kg	50 kg

Specifications	Diana II ELA, 90556 with Hanger Bar 2P	Diana II ELA, 90556 with Hanger Bar 4MSB
Maximum lifting capacity	160 kg	160 kg
Rotation of the hanger bar	360 degrees	360 degrees
Rechargeable battery (incl. charging cable)	Capacity 2.9 Ah	Capacity 2.9 Ah
Frame material	Powder coated steel	Powder coated steel
Emergency lowering function	Yes	Yes
Hand control functions	Up/down, leg spreading, battery status, service indicator	Up/down, leg spreading, battery status, service indicator
Certification	CE, ISO 10535	CE, ISO 10535



*We improve freedom of movement
and quality of care for people.*

www.humancaregroup.com



HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
+46 8 665 35 00
info@humancaregroup.com

HUMAN CARE GERMANY

Industriestraße 13, Gebäude 1/9
63755 Alzenau
+49 (0) 6023 994 936 8
info.de@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
ON K2E 7K1 Ottawa
+1 613 723 6734
info.ca@humancaregroup.com

HUMAN CARE THE NETHERLANDS

Elspeterweg 124
8076 PA Vierhouten
+31 577 41 21 71
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron rd, Suite K
TX 78754 Austin
+1 512 76 7199
info.us@humancaregroup.com