



Reliable  
Maximum user weight  
230 kg  
Large range



## For comfortable and safe transfers

The Diana II Comfort is a stable, electrically operated lift with a high maximum user weight, allowing people weighing up to 230 kg to be safely supported when they are being lifted. Our extensive range of E-clip slings provide the required support.



**Lifting from a bed**, in/out of a chair or wheelchair, on/off the toilet and even from the floor, these are a few examples of the many situations where the Diana II Comfort can provide support. This lift is suitable for clients of different heights and weight, due to the high load capacity, the large range of the lift arm and the adjustable 4-point spreader bar.

**The electric leg operation** allows the lift to be placed around a chair/toilet to get as close as possible to a client for a more comfortable transfer. It is easy to maneuver the lift thanks to the large, high-quality twin casters. This, and many other factors, make the Diana Comfort II, a durable high-quality product.



Using the 4-point spreader bar allows you to place the client into the appropriate recumbent or sitting position. A practical solution is the manually operated spreader bar: it can provide three positions: recumbent, semi-recumbent and sitting. For frequent use, we recommend an electric 4-point hanger bar, where every angle inbetween recumbent and sitting is a possibility. It has a rotation point so that the client can easily and safely be rotated 360 degrees. This greatly increases the level of user convenience.

**Optimal support** can be provided by combining the Diana Comfort II with one of our slings. Our slings are attached using the unique E-clips. The universal E-clip is a proven system for fastening a sling to a lift, thereby significantly reducing the risk of errors. Our slings provide optimal support for extra secure, comfortable and safe transfers.



Mobility Aids



Healthcare Beds



Lifting Solutions



Hygiene



Convertible Chairs

## Accessories



**Full Body Sling E-clip Soft**  
Provides full support and helps to achieve a comfortable and safe transfer of the client. The sling is attached to the lift by means of our E-clips.  
**25010C-X (25005C-X in Net)**



**Columbus Sling E-clip**  
The same as the Full Body Sling and even provides support under the client's sitting area.  
**Article number 25040C-X**



**Toileting Sling E-clip**  
Comfortable and sturdy toileting sling. Ideal for transporting clients to and from the toilet.  
**Article number 25085C-X**



**Soft padding for the hanger bar**  
Option to decide where to place the padding on the hanger bar.  
**Article number 50866**

## Use



**Lifting in/out of bed**  
Provides support to clients during lifting from the bed and can help with positioning on the bed.



**From and to a chair**  
Suitable for transferring to and from a chair or wheelchair and positioning the client in a chair.



**Lifting from the floor**  
Provides support to clients while being lifted from the floor.



**Toileting**  
Suitable for transfer to and from the toilet.

## Features



**Casters**  
2 twin castors, with brakes (Ø 125 mm, rear) and 2 twin castors, without brakes (Ø 100 mm, front).



**Hand control**  
User-friendly handset, easy to use, ergonomic design with durable coiled cable.



**Locking**  
Lock for the rotation of the electrically operated 4-point spreader bar for extra safety.



**Battery**  
External rechargeable battery with Linak components.

## Technical specifications

Specifications	Diana II Comfort 4ESB 90560	Diana II Comfort 4MSB 90559
Length of the chassis	1300 mm	1300 mm
Width of the chassis	670 mm	670 mm
Height of the mast	1350 mm	1350 mm
Min. Width leg spreading (internal/external)	520 – 670 mm	520 – 670 mm
Max. Width leg spreading (internal/external)	1030 – 1180 mm	1030 – 1180 mm
Height of the push handle	950 – 1270 mm	950 – 1270 mm
Turning radius	1440 mm	1440 mm
Diameter of the casters	125 mm (rear, with brake), 100 mm (front)	125 mm (rear, with brake), 100 mm (front)
Height of the frame	115 mm	115 mm
Lifting range	680 – 1960 mm	680 – 1960 mm

Specifications	Diana II Comfort 4ESB 90560	Diana II Comfort 4MSB 90559
Max. load capacity	230 kg	230 kg
Product weight	66 kg (incl. 4ESB)	66 kg (incl. 4MSB)
Tilting angle of the hanger bar	-17 to + 42 degrees	13, 29, 45 degrees
Rotation of the spreader bar	360 degrees with locking	360 degrees
Rechargeable battery (incl. charging cable)	Capacity 2.9 Ah	Capacity 2.9 Ah
Frame material	Powder coated steel	Powder coated steel
Emergency lowering function	Yes	Yes
Hand control	up/down, leg operation, battery status, service indicator	up/down, leg operation, battery status, service indicator
Certification	CE, ISO 10535	CE, ISO 10535



*We improve freedom of movement and quality of care for people.*

[www.humancaregroup.com](http://www.humancaregroup.com)

