



Ergonomic
Max. user weight 230 kg
Compact



Perfect for active transfers

The seBastian II is a compact, electrically operated stand-up lift with a high maximum user weight, allowing people weighing up to 230 kg to be safely supported while getting up or standing up. Extra support and comfort can be provided by using our Back Belt E-clip.



Multifunctional and durable

This fully electric lift provides the client with support whilst carrying out daily activities, such as getting out of bed, sitting down in a chair, going to the toilet, dressing and undressing, as well as during rehabilitation exercises.

In addition to its high load capacity (up to 230 kg), the seBastian II also has a high lifting range. This lifting range together with the height-adjustable lower leg supports means that this lift is suitable for clients of different heights.

makes it easier to get closer to a client and therefore achieve a more comfortable transfer.

Optimal support

The seBastian II can be fitted with our unique Back Belt E-clip sling. The universal E-clip is a proven system for fastening a sling to a lift and thereby significantly reducing the risk of human error.

Our sling provides optimal support for extra secure, comfortable and safe transfers. The sling encourages a natural movement because the padded back provides a comfortable support in the client's back during the upward direction of movement. The Back Belt E-clip also ensures the client has a secure feeling and prevents any risk of falling backwards. Use of the Back Belt E-clip means that the caregiver has his/her hands free during the transfer process.



This extremely compact lift is very easy to move and maneuver thanks to its compact size and large size high-quality twin castors.

Due to the electrically operated extra wide leg spreading, the seBastian II can easily be placed around a chair or toilet. This



Accessories



Back Belt E-clip
Provides back support for greater comfort and safety. Attaches using secure E-clips.
Article number 25130C-(X)



Battery
Extra Linak battery to ensure the lift is always ready for use.
Article number 04-N0101-0012



External charger
Charger for external battery charging.
Article number 04-N0105-0018

Use



Lifting in/out of bed
Provides support to clients during lifting from the bed and can help with positioning on the bed.



From and to a chair
Suitable for transferring to and from a chair or wheelchair and positioning the client in a chair.



Toileting
Suitable for transfer to and from the toilet.

Features



Casters
2 twin castors with brakes (Ø 125 mm, rear) and 2 twin castors without brakes (Ø 100 mm, front).



Hand control
User-friendly handset, easy to use, ergonomic design with durable coiled cable.



Footplate
Horizontal anti-slip footplate to keep the feet in place while standing up/sitting down.



Lower leg supports
Height-adjustable, bowl-shaped lower leg supports with built-in angle adjustment that follows the standing movement for extra comfort.

Technical specifications

Specifications	seBastian II, 90554
Length of the chassis	1080 mm
Width of the chassis	658 mm
Width of the lift (including push handles)	750 mm
Height of the mast	1130 mm
Minimum width of the chassis (internal - external)	510 – 660 mm
Maximum width of the chassis (internal - external)	1010 – 1150 mm
Lifting range	780 – 1740 mm
Turning radius	1240 mm
Height of the chassis	115 mm

Specifications	seBastian II, 90554
Diameter of the casters	125 mm (rear, with brake), 100 mm (front)
Total weight	60 kg
Maximum lifting capacity	230 kg
Rechargeable battery (incl. charging cable)	Capacity 2.9 Ah
Frame material	Powder coated steel
Emergency stop	Yes
Emergency lowering function	Yes
Hand control functions	Up/down, leg spreading, battery status, service indicator
Certification	CE, ISO 10535



We improve freedom of movement and quality of care for people.

www.humancaregroup.com

