


# Alpha Turning Disc

Durable  
Easy turning  
Visible



The Alpha Turning Disc can make transfers between bed and (wheel)chair easier for clients that are able to stand. The start and stop function of the Alpha Turning Disc allows movements to be precisely controlled.

We recommend the use of a Turning Disc in combination with the Alpha Footstool.

## Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Turning Disc	PM-5080	40 cm diameter / 15,7" diameter	ABS coated plastic fabric, latex free	200 kg / 440 lb



CE MD = 200 kg / 440lb

**PETERMANN**  
part of the Human Care Group

These products are distributed by  
Human Care AB for Petermann GmbH.



Visit [www.humancaregroup.com](http://www.humancaregroup.com) for the latest product information and contact details to national offices and partners.