

Flexible Sliding Board

Smooth
Flexible
Sideward
transfer



PM-5050 / PM-5020



PM-5030



The Alpha Flexible Sliding Board has a very smooth surface in order to minimize resistance due to friction. It is used to transfer clients, e. g. from the edge of a bed to a wheelchair, toilet or car seat. The client does not have to be able to stand up, but can make the change in position while sitting, either with or without assistance.



Technical specification

Product Name	Art. nr.	Dimensions	Material
Alpha Flexible Sliding Board Small	PM-5020	58 x 32 cm / 22,8" x 12,6"	Polyethylene, latex free
Alpha Flexible Sliding Board Medium	PM-5050	74 x 32 cm / 29,1" x 12,6"	Polyeten, latexfree
Alpha Flexible Sliding Board with rounded corners	PM-5030	74 x 32 cm / 29,1" x 12,6"	Polyethylene, latex free



 = 200 kg / 440lb

PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.