

Alpha Non Slip Mat

Prevents sliding in bed



The Alpha Non Slip Mat prevents the feet from slipping in bed. When placed underneath the client's feet in bed, the non-slip mat helps to prevent slipping while repositioning or raising the patient.

A non-slip mat increases the frictional resistance by 100 %. It also makes it possible for the client moving toward the top of the bed by pressing with their feet.

Alpha Non-Slip Mat on Roll is a non-slip plastic film with high friction that comes on a roll. Cut the section or shape you want and let the product meet several anti-slip needs.



Technical specification

Product Name	Art. nr.	Dimensions	Material
Alpha Non Slip Mat Small	PM-7010	35 x 25 cm 13,7" x 9,8"	Polymeric Compound, red
Alpha Non Slip Mat Medium	PM-7011	45 x 38 cm 17,7" x 14,9"	Polymeric Compound, red
Alpha Non Slip Mat on Roll	PM-7007	40 x 200 cm 15,7" x 78,7"	Polymeric Compound

