



Alpha Lifting Seat

Safe
Reliable
Stable



The Alpha Lifting Seat is ideal for stable seated positioning when transferring persons down staircases and other transfer situations where there is no room for mechanical aids. The weight of the client can be distributed over up to five persons, even of different heights, to save their backs. The stable safety handles make the client transfer safe and easy.



Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Lifting Seat	PM-4035	38 x 43 x 55 cm 14,9" x 16,9" x 21,6"	Nylon	180kg / 396lb




 = 180 kg / 396lb

PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.