


# Alpha Stand-up Aid Quick

Universal  
Many applications  
Supportive



The Alpha Stand-up Aid Quick is a universal aid with many potential applications. It is a safe transfer aid, which allows for a transfer of the client in the bed or from the bed to the wheelchair. It supports the client in a pleasant and comfortable manner.

The Alpha Stand-up Aid Quick provides support on the client's back. This application invites the client to lean forward in a physiological manner to facilitate a natural motion.

This product has 4 handles on both sides in order to be easy-to-use and back-friendly for the caregiver. It can also be used in combination with our Flexible Sliding Board.

## Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Stand-up Aid Quick	PM-6080	145 x 11 cm / 57" x 4,3"	Polyester	200 kg / 440 lb



**PETERMANN**  
part of the Human Care Group

These products are distributed by  
Human Care AB for Petermann GmbH.