

Alpha Turning Sheet & Alpha Fitted Sheet

For pain sensitive, heavy or special care clients




The Alpha Turning Sheet and the Alpha Fitted Sheet work great together in combination: a large fitted sheet which is placed on the mattress and a large sliding sheet. The upper side of the sliding sheet is made of cotton/natural ticking.

The Turning Sheet is used for pain-sensitive clients, heavy clients, and special care clients. Moving and laterally repositioning clients is made much easier. Rolling up the Turning Sheet against the clients body in a lateral position and then fixing it in place with a bed sheet provides simple and gentle lateral position support.



Technical specification

Product Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Fitted Sheet	PM-1584	90 x 200 cm 35,4" x 78,7"	Cotton, latex free	200 kg 440 lb
Alpha Turning Sheet	PM-1582	140 x 200 cm 55,1" x 78,7"	Cotton, natural thicking, latex free	200 kg 440 lb



These products are distributed by Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.