


Alpha Turning Sheet with handgrips

Breathable
Soft
Hidden
handgrips



This Alpha Turning Sheet with handgrips was developed specially for clients who are particularly sensitive to pain. The processing of breathable textile products enables the client to remain on the Turning Sheet on an ongoing basis.

The hidden handles and loops of the turn sheet make it possible for the carer to move the client freely and gently to the desired care situation in the bed. The Turning Sheet can be attached to the head of the bed or laterally to the slatted frame with the attachment straps provided.

Technical specification

Article Name	Art. nr.	Dimensions	Material	Max lifting weight
Alpha Turning Sheet with handgrips	PM-2098	100 x 200 cm 39,3" x 78,7"	Cotton, latex free	200 kg 440 lb



CE MD = 200 kg / 440lb

PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.