HUMAN CARE
Your Life. Your Way.




Sliding Sheets

- for easier positioning, gentle care and reduced strain.

Sliding and drawing sheets from Human Care

Positioning in bed is important to increase comfort and reduce the risk of pain and tissue damage. At the same time, repositioning can be challenging for the user and physically demanding for caregivers. Our sliding and drawing sheets reduce friction between the patient and the surface, making moving and positioning in bed smoother and more comfortable. They are available in several sizes and models for different needs.



Sliding sheet - Alpha Base Sheet 2-way


The Alpha Base Sheet 2-way sliding sheet has a low-friction satin surface that makes it easier to move and turn around the bed. It is available in several sizes and models. For patients with greater mobility, this sheet may be enough to easily turn and position themselves in bed.

- Shorter versions are available for more active users
- Full-length models are available for less mobile users
- A form-fitting version is available that also fits airborne antidecubitus mattresses

There are also sliding sheets where the sliding surface has an incontinence function - Alpha Base Sheet Waterproof.

Drawing sheet - Alpha Draw Sheet

For patients who need more help positioning themselves in bed, the sliding sheet can be combined with a drawing sheet. The drawing sheet is placed on top of the fitted sheet and is available in two designs.





Drawing sheet - Alpha Draw Sheet 2-way

This is the simpler variant of the draw sheets that are placed on top of the sliding sheet with the striped lines placed horizontally across the width of the bed, from side to side. It is intended for lateral movement and positioning/turning and helps when assistance from caregivers is required. Caregivers can then help by pulling on the sheet.







Sliding Draw Sheet - Alpha Draw Sheet 4-way

Allows for movement and positioning with minimal friction both sideways and also up and down in bed. This makes the Alpha Draw Sheet 4-way particularly suitable for heavier users or users with reduced mobility and need extra support when repositioning. On the back, there is a sliding fabric that provides low friction, which increases the possibility of efficient multi-directional movement. The crossed lines on the sheet show that it has a sliding function in four directions.

Available with integrated handles that can be used manually by caregivers or connected to a hanger bar for additional relief.

The model with handles is also available in incontinence version for larger amounts of fluid - Alpha Draw Sheet 4-way Handles Inco.



- **Low friction – easy positioning with minimal force**
- **Easy to handle – intuitive for caregivers, empowering for the user**
- **High comfort – soft glide surface reduces pressure and friction**
- **Modern design – a stylish form that blends into any environment, in harmonious colours.**

Sliding sheet - Alpha Base Sheet 2-way

Article number	Product name	Dimensions – product (length × width)	Dimensions – Sliding surface (length × width)	Fits bed width
PM-1611	Alpha Base Sheet 2 way 100×200 cm	100×200 cm	100×50 cm	80 cm
PM-1629	Alpha Base Sheet 2 way 100×200 cm	100×200 cm	100×60 cm	80–90 cm
PM-1616	Alpha Base Sheet 2 way 100×220 cm	100×220 cm	100×75 cm	105 cm
PM-1667	Alpha Base Sheet 2 way 100×220 cm	100×220 cm	100×220 cm	120 cm
PM-1612	Alpha Base Sheet 2 way 140×200 cm	140×200 cm	140×50 cm	80 cm
PM-1634	Alpha Base Sheet 2 way 140×200 cm	140×200 cm	140×60 cm	80–90 cm
PM-1621	Alpha Base Sheet 2 way 140×220 cm	140×220 cm	140×75 cm	105 cm
PM-1668	Alpha Base Sheet 2 way 140×220 cm	140×220 cm	140×95 cm	120 cm
PM-1613	Alpha Base Sheet 2 way 200×200 cm	200×200 cm	200×50 cm	80 cm
PM-1639	Alpha Base Sheet 2 way 200×200 cm	200×200 cm	200×60 cm	80–90 cm
PM-1622	Alpha Base Sheet 2 way 200×220 cm	200×220 cm	200×75 cm	105 cm
PM-1669	Alpha Base Sheet 2 way 200×220 cm	200×220 cm	200×95 cm	120 cm
PM-1681	Alpha Base Sheet 2 way 200×280 cm	200×280 cm	200×60 cm	140–200 cm
PM-1641	Alpha Base Sheet 2 way fitted 200×90 cm	200×90 cm	200×70 cm	80–90 cm
PM-1651	Alpha Base Sheet 2 way fitted 200×105 cm	200×105 cm	200×80 cm	105 cm
PM-1661	Alpha Base Sheet 2 way fitted 200×120 cm	200×120 cm	200×95 cm	120 cm
PM-1629W	Alpha Base Sheet 2 way waterproof 100×200 cm	100×200 cm	100×60 cm	80–90 cm
PM-1634W	Alpha Base Sheet 2 way waterproof 140×200 cm	140×200 cm	100×60 cm	80–90 cm

Material: Cotton/terylene (C/T 55/45%), Cotton/terylene (C/T 55/45%) and terylene, polyuretan

Drawsheet - Alpha Draw Sheet 2-way


Article number	Product name	Dimensions – product (l×b)	Recommended bed width
PM-1631	Alpha Draw Sheet 2 way 100×200 cm	100×200 cm	80–120 cm
PM-1636	Alpha Draw Sheet 2 way 140×200 cm	140×200 cm	80–120 cm

Material: Cotton/terylene (C/T 55/45%)

Sliding Draw Sheet - Alpha Draw Sheet 4-way

Article number	Product name	Dimensions – product (l×b)	Recommended bed width
PM-1637	Alpha Draw Sheet 4 way 140×200 cm	140×200 cm	80–120 cm
PM-1647	Alpha Draw Sheet 4 way 200×200 cm	200×200 cm	80–120 cm
PM-1657	Alpha Draw Sheet 4 way Handles 140×200 cm	140×200 cm	80–120 cm
PM-1687	Alpha Draw Sheet 4 way Handles 200×200 cm	200×200 cm	80–120 cm
PM-1649	Alpha Draw Sheet 4 way Handles Inco 140×120 cm	140×120 cm	80–120 cm

Terylene polyester, viscose, polyurethane layer, terylene, polypropylene handle



Ironing at low temperature is allowed

<200 kg

CE

MD

CERTIFIED ISO 13485 MEDICAL DEVICES • QUALITY MANAGEMENT

HUMAN CARE SWEDEN
Årstaängsvägen 21B, 117 60 Stockholm
Tel: +46 8 665 35 00
Email: info.se@humancaregroup.com

98944, EN, Rev1, Alpha slide sheet, Prd sheet, Aug26